DCC Chromatograph B-670

97033 500 382

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Installation

Do not place the B-670 DCC chromatograph in direct sunlight, in order to avoid a shift of the phase equilibrium and expansion of the solvent system.

Separation of temperature sensitive samples may be done in a refrigerated room. However, condensation may occur. The instrument should not be kept in a refrigerated room for a long time to minimize condensation and eventual corrosion.

Before you turn the instrument on for the first time

Check that voltage and frequency are the same as the marks on the type plate.

For **final inspection** and to fulfill shipping regulations, the columns have been filled with **Tetrachlorethylene**.

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Selection of solvents

Reprinted with the permission of the author:

Prof. Dr. K. Hostettmann

A quick way to select a solvent system consists in checking the sample by thinlayer chromatography (TLC) on silica gel with the water-saturated organic layer [11, 12]. Empirically, we found that if the R_f values of the compounds to be separated are higher than about 0.50 (less polar solutes), the less polar layer is suitable for use as the mobile phase. With more polar substrates ($R_f < 0.50$), the more polar layer should be used as the mobile phase. If the R_f values are in the range 0.40–0.60, the separation can be achieved by using either the more polar layer or the less polar layer as the mobile phase. This is illustrated in the following by the separation of a fraction (60mg) of a crude extract of a Rocky Mountains Gentiana species (*Gentiana strictiflora* [RYDB.] A NHS) obtained after column chromatography on Sephadex LH-20. On silica gel TLC, with $\text{CHCl}_3\text{-MeOH-PrOH-H}_2\text{O}$ (9:12:1:8), lower layer, this fraction showed two UV active spots at $R_f = 0.35$ and 0.41. DCCC separation by using the more polar upper layer as mobile phase yielded two pure compounds 1 and 2. The use of the more polar layer as the mobile phase resulted in earlier elution of the more polar compound 1. For comparison, the same mixture of 1 and 2 was submitted to DCCC and eluted with the less polar lower layer. The resolution was better than in the previous run. As expected, the less polar glycoside 2 was eluted first. These comparative experiments suggest that in the case of difficult separations, it could be better to use the less polar layer as mobile phase for eluting compounds with relatively low R_f values, and the more polar layer for eluting compounds with relatively high R_f values. A better resolution will be obtained, but the separation time will be much longer. Solvent system selection based on TIC is very simple and rapid and gives generally satisfactory results. But since DCCC depends solely on the partition coefficient, separation cannot always be predicted from the TLC behavior.

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

[1] K. Hostettmann, *Planta med.* 39, 1 (1980)

DCCC separation of xanthone-O-glycosides with $\text{CHCl}_3\text{-MeOH-n-PrOH-H}_2\text{O}$ (9:12:1:8). Top elution curve: the mobile phase was the more polar upper layer. Bottom elution curve: the mobile phase was the less polar lower layer.

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- ① Power On/Off
- ② Pump On/Off
- ③ Pump Motor Control
- ④ Pressure Warning
- ⑤ max. Pressure set

- ⑥ Sample Chamber
- ⑦ Ascending/Descending
- ⑧ Vent Valve
Vent Outlet
- ⑨ Pump Stroke

Controls: Front panel

Description

- ① Power on/off
Main switch with built-in bulb.
Power supply for the pressure sensing amplifier print GS 2.
Resets the electronic lock of the pressure monitor.
- ② Pump on/off
Switch with built-in bulb.
Power supply for the print GS 1, motor control
- ③ Pump Motor Control
Potentiometer for the speed regulation of the pump motor.
- ④ Pressure Warning
This bulb is on, whenever the pump is switched off by the automatic pressure monitor.
By pressing the mains switch (Power on/off), the electronic lock will be reset.
- ⑤ max. Pressure set
Potentiometer to set the pressure limit to shut off the pump.

- ⑥ Sample Chamber
The volume of the sample chamber is 15 ml. The inlet end is marked with an arrow. To avoid overloading of the first columns, the sample is initially pumped into the sample chamber.

The sample solution (1 : 1 mixture of mobile and stationary phase) is collected in the sample chamber. There the sample solution is slightly diluted with stationary phase and a new equilibrium under separation conditions will be established. The Sample chamber can be rotated so that during a separation only the mobile phase is in the sample chamber and a change of phase can be made in the sample chamber.

- ⑦ Ascending/Descending
With these connections and according to the flow diagram, the desired mode of operation is established.

- ⑧ Vent Valve/Vent outlet
To purge air from the system and to facilitate the suction when starting the pump.

- ⑨ Pump Stroke
Mechanical control of the plunger stroke length of the metering pump.

Caution

Never set stroke beyond O mark. Adjusted beyond zero, plunger may be damaged.

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Tab. I

Filling of columns with stationary phase

The columns should be filled during the night so that the instrument is ready for a separation on the following morning.

Read the required volume of stationary phase from table I and fill the corresponding amount into the flask ⑳.

When necessary replace the filtering material (glass wool) of filter tube ⑯.

Insert aspirating tip ⑩ into the flask ㉑.

Mains ① «on».

Potentiometer ⑤ on 40bar.

Check that the connections ⑦ (ascending/descending) are properly connected for the desired operation mode.

Pump ② «on».

If changing the stationary phase rotate the sample chamber to displace the previous phase.

The pump might not pull any of the phase if air is in the system. Open vent valve ⑧ to purge the air, or aspirate with a syringe (to be connected at «Vent outlet») phase through the open valve.

The time to fill the columns can be estimated by reading the volume (from table I) and the flow rate (from nomogram I).

Adjust chosen flow rate with potentiometer ③ and knob ④.

Caution

Never set knob ④ beyond the O-mark. Adjusted beyond the zero mark, the plunger could «bottom» the pump head and break.

Fit the outlet ⑮ to a flask, fraction collector or detector.

Set the pressure, at which the pump will be turned off by the pressure monitor, with potentiometer ⑤.

Columns can be viewed by opening the swivel door ㉒ and pulling out the first rack ㉓. The bar ㉔ must be secured to prevent a complete, accidental withdrawing of the rack.

Sample conditioning

The samples must not contain any suspended matter and/or particles. Filtering of any type may be necessary.

Note: Clogged capillary tubes will result if samples contain any foreign particles.

Separations of crude extracts will contaminate the system quicker than separations with semi purified samples.

Volumes of the various racks and columns

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
Loading the sample

Theory requires a punctual sample introduction at the beginning of the first column and is therefore enhanced. Practically this can not be achieved due to the sample's solubility. The volume of the sample solution must be kept as low as possible. However a to a high concentration may change the equilibrium of the phase system. Depending on the separation problem (solvent system, sample quantity, volume of sample solution), the use of a smaller sample chamber may be indicated. The standard sample chamber has a volume of 5 ml. The smaller ones will reduce dilution of the sample solution. Also the possible resolution will be increased (less tailing/bleeding).

1 ml sample chamber 16545
15 ml sample chamber 16295

- Dissolve the sample into a 1 : 1 mixture of mobile and stationary phase. Place the sample tube into the holder ④.
- Check position of sample chamber ⑥. The arrow must point in the opposite direction of the desired flow of droplets.
- Pump ② «off».
- Dip aspirating tip ⑩ into sample. If necessary use filter support ⑧.
- Pump ② «on».
- Pump ② «off» before air is aspirated.
- Reduce flow rate (potentiometer ③ and knob ⑨). 15 to 50 ml/h are typical flow rates of the mobile phase. Depending upon the separation problem, the column diameter, the viscosity of the mobile phase a.s.f., the flow rate may be higher or lower than indicated.
- Dip aspirating tip ⑩ into mobile phase. Mix the sample again by rotating sample chamber.
- Pump ② «on».
- Observe and optimize generation of droplets with potentiometer ③ and knob ⑨.

Note!

Too large of a drop of mobile phase will displace the stationary phase. This is avoided by keeping the flow rate reasonably low.

Many small droplets will give better separations than a few large drops!

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
UV-Detector

Observe the following points if a UV-detector is used together with the B-670 DCC chromatograph:

- Cell

As a general recommendation the flow through cell should have a small optical pathlength (0.5 or 1 mm for 20 to 50 ml/min) since the sample quantities in DCC Chromatography are high. Also the detector should have a time constant of approx. 10 seconds to obtain a neat recording. Such a time constant may be made as illustrated:

- Accumulation of stationary phase in the cell must be avoided. Depending upon the operation mode (ascending/descending), rotate the detector 180°C.
- Bubble formation can occur because expansion causes gas bubbles to form. Increase the pressure in the cell at the outlet by adding a length of teflon tubing (part 16 543).

Refractive index detector

The use of this type of detector is not possible with your DCC chromatograph since a minimum amount of stationary phase is always displaced with the mobile phase.

Change of phase

After a separation is completed aspirate 20 ml methanol as <buffer> and than the remaining stationary phase is displaced with the new one needed for the next separation. Rotate the sample chamber to displace the previous mobile phase.

If the instrument will not be used for another separation, Methanol is filled up to the vent valve (⊕) to prevent sticking of the ball valves to the seats.

Cleaning the columns

It is necessary to clean the columns when no more droplets are generated due to contamination; the droplets <adhere> to the column wall.

The following steps must be strictly followed:

1. Clean filter (Ⓜ).
2. Aspirate 20ml Methanol.
3. Aspirate 50ml H₂O dist.
4. Aspirate over night a 20% Extran solution (Merck 7555).
5. Rinse with distilled H₂O to pH7.
6. Aspirate 20ml Methanol.
7. Fill in new stationary phase.

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
Changing columns and racks

1. Withdraw securing bar ③①.
2. Swivel out rack holders ③②.
3. Hang rack on holders.
4. Remove springs ③③ of column securing plate ③④.
5. Take column off at the couplings and remove it in upwards direction.
6. Remove teflon tubes from their guides and place column on the bench.
7. Open coupling with wrenches.
8. The nut ③⑦ is hold on the column flange by means of the plastic «spaghetti» ③⑧. A small pair of tongs are well suited to extract and to insert the «spaghetti» ③⑧.

How to make teflon flanges


1. Cut the teflon tube with a razor blade as shown.
2. Fit the components of the coupling onto the teflon tube.
3. Grip the teflon tube with the flanging tool ④⑩ be sure that 1–2mm of the tube are projecting the tool.
4. Press the hot flange soldering iron ③⑨ firmly against the tube.

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




Drawer liquid section

⑥ Sample chamber 15 ml	16 295
1 ml sample chamber	16 545
5 ml sample chamber (Standard)	16 546
⑦ Coupling block 3F	16 352
⑧ Vent valve	16 262
⑩ Aspirating tip	16 342
⑪ Coupling	16 296
⑫ Coiled tube	16 341
⑬ Metering pump	16 331
⑭ Pressure sensor	
⑮ Outlet coupling 1F	16 353
⑯ Filter tube	16 491
⑳ Window	16 323
㉑ Bottom cover	11 336
㉒ Sample holder	16 360

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Components and Spare parts, Drawer liquid section

㉓ Pear flask 1l	16 326
㉔ Stopper	16 325
⑩ Aspirating tip	16 342
⑫ Coiled tube	16 341
⑦ Coupling block 3F	16 352
⑭ (page 29) Pressure sensor	
⑧ Vent valve	16 262
Needle for 16 262	16 257
Gasket for 16 262	16 259
⑥ Sample chamber	16 295
⑬ (page 23) Metering Pump	
⑯ Filter tube	16 491
㉕ (page 25) Motor	
Teflon tube 5 m	16 543

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


Coupling for: ⑥ **Sample chamber,** ⑩ **Aspirating tip,**
 ⑬ **Filter tube,**
 ⑥② Nut M 9,5 spec. 16 296
 ⑥③ Spaghetti 16 324

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Coupling for: ⑩ **Aspirating tip,** ⑫ **Coiled tube,**
 ⑬ **Pump inlet,**
 ⑥⑨ Coupling complet 16 544

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Vent valve ⑧

1. Loose fastening nut on front panel and open the valve cover by unscrewing the three countersink bolts.
The gasket may now be changed.
2. When assembling the vent valve, the valve needle must be removed from the valve cover.

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Sample chamber ⑥

1. Unscrew the lateral countersink bolts.
2. After loosening the «allen» screws, the holders can be taken off from the coupling.
3. The tube coming from the pump, must be connected to the end of the sample chamber which is marked with an arrow.

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Description of the metering pump

Basically the pump consists of the drive unit (17), the plunger (18) and the pump head (19).

The flow rate may be varied either with the knob (9) which changes the stroke length, or by varying the speed of the drive unit (17).

Plunger (18) and ball valves (20) are made of sapphire.

Care

Never set stroke (9) beyond 0-mark. Adjusted beyond zero, plunger may be damaged.

Maintenance and care of the metering pump

1. Plunger gasket (21)

1.1 Tightening of the gasket is done by tightening the nut (22) with the rod, supplied with the instrument.

1.2 Replacing the gasket (21):

1. Set knob (9) on 100 to retract the plunger (18) to its rearmost position. This decreases the risk of the plunger breaking.
2. Loosen the knurled nuts (23) with alternating ½ turns.
3. Carefully remove the pump head (19).
4. Open the nut (22) with the rod supplied with the instrument and remove the gasket (21) consisting of 5 parts.
5. Insert new gasket set observing orientation and order of the gasket parts strictly.
6. Set knob (9) on zero.
7. Assembling is done in the reverse order. Dust and particles will decrease the life of the new gasket and plunger.
8. Observe the flow direction when mounting the pump head.

2. Replacement of ruby plunger (18).

When changing the plunger, the gasket (21) must be replaced too.

- 1.–6. Refer to replacement of gasket (21).
7. Withdraw plunger (18) and open guiding screw (24) to take away the plunger. Do not scratch surface of metal piston (42).
8. Mount new plunger.
9. Clean and oil metal piston (42) and bearing.
10. Set knob (9) on zero and mount assembly in reverse order.

3. Valves (43) and (44) of metering pump

The valves (43) and (44) consist of identical parts. The difference is the direction of mounting of the double ball valves (25).

When assembling the valves (43) and (44), the flow direction must be observed.

4. Replacement of double ball valve assembly

1. When changing the double ball valve assembly, gasket (26) must be replaced too.
2. Unscrew valve cartridge (27) and remove all components.
3. Observe flow direction during assembling! Dust and particles will lead to leakage of the ball valves.

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
Print GS 2 und 3

**Diagram print GS 2, Pressure measuring amplifier
PA Pressure sensor
GS 3 Gauge adjustment
Gauge**

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
For calibration reasons, the pressure measuring amplifier and pressure sensor are only available as one part. This also applies for the instrument and its adjustment print.

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
	New	Standard Exchange
GS 2 Pressure measuring amplifier and PA pressure sensor	13 667	16 535
Gauge adjustment print with gauge	13 668	-

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Components and Spare parts, electronics

⑧1	Gauge with adjustment print	13 668
⑧2	Potentiometer 5 kΩ	13 247
	Knob	13 759
	Cover	13 761
⑧3	Indicating ring	13 782
	Lamp	01 594
	Cover, red	01 612
⑧4	Bulb 5 pcs.	01 559
	Press button	01 582
	Cover, yellow	01 597
⑧5	Bulb 5 pcs.	01 559
	Press button	01 288
	Cover green	01 614
⑧6	Bulb 5 pcs.	01 559
	Condenser 5 μ F	13 181
⑧7	Transformer	13 688

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Trouble shooting

Trouble	Possible cause	Remedy
Lamp «on/off» does not light.	<ol style="list-style-type: none"> 1. Main fuse blown. 2. Power cord. 3. Fuse F1 blown. 4. Fuses F 2/F3 blown. 5. Bulb faulty. 6. Switch ① faulty. 	<ol style="list-style-type: none"> 1. Check and replace. 2. Check and plug in power cord. 3. Replace fuse at rear panel, check for correct voltage. 4. Check and replace fuses in electronics compartment. 5. Replace bulb. 6. Replace switch.
Lamp «Pressure warning» is light.	<ol style="list-style-type: none"> 1. Pressure monitoring circuitry has turned the pump off, due to too high pressure. 2. System clogged. 	<ol style="list-style-type: none"> 1. Set Potentiometer ⑤ at a higher value. Reset electronic lock by pressing switch ① off and on. If the pump is repeatedly being switched off, reduce the flow rate. 2. Particles or dirt have clogged filters, connectors or columns. To locate the obstruction quickly, pressure is maintained on the system (switch ① on) By observing the pressure gauge, open connections from the outlet forward, opposite to the flow direction. A sudden drop of the pressure means that the obstruction is between this connection and the one opened before. To locate the obstruction, the system has to be pressurized repeatedly until the pressure monitoring circuitry switches the pump off.
Pump motor does not run.	<ol style="list-style-type: none"> 1. Fuse F2/F3 blown. 2. Print GS 1, GS 2 or motor faulty. 	<ol style="list-style-type: none"> 1. Check and replace fuses in the electronics compartment. 2. Check whether relays A and B on print GS 2 are switching. If they operate correctly, the print GS 1 or the motor are faulty.
Pressure indication not working	<ol style="list-style-type: none"> 1. Sensor PA or print GS 2 faulty. 2. Gauge faulty. 	<ol style="list-style-type: none"> 1. For calibration reasons, these parts have to be replaced together. 2. Replace gauge and its adjustment print.
Phase not being pumped	<ol style="list-style-type: none"> 1. Air in the system 2. Plunger sticking, is not pushed backwards. 3. Leak. 4. Ball valves sticking. 5. Ball valves dirty. 6. Plunger gasket leaking. 	<ol style="list-style-type: none"> 1. Open vent valve to purge and to facilitate suction action of pump, or aspirate phase with a syringe (to be connected at «Vent outlet») through the vent valve. If air is already in the columns, change connection ascending/descending to displace the air. 2. Push plunger at the guiding screw ⑳ backwards. Motor must be switched on. 3. Check all connections, starting at the aspirating tip, for leakage. 4. Loosen the ball valves with gentle knocks against the ball valve body. 5. Remove ball valve assembly and clean it from particles and other dirt. 6. Tighten or replace gasket.
No droplet generation with proven phase system	<ol style="list-style-type: none"> 1. Columns contaminated. 2. Pressure too high. 3. Phase system. 4. Operation mode. 5. Wrong phase. 6. Wrong columns. 7. Equilibrium of phase system. 	<ol style="list-style-type: none"> 1. Clean columns exactly according instructions. 2. Reduce flow rate and therefore the pressure. 3. Make sure that the phase system has reached equilibrium (mix it well and allow more time). 4. Verify that the desired operation mode is possible with this phase system. 5. Check if the phases have been mixed up. 6. Droplet generation is possible with different diameter columns. 7. A too high sample quantity has shifted the equilibrium.
No droplet generation with new phase system	<ol style="list-style-type: none"> 1. As above. 2. Column diameter. 	<ol style="list-style-type: none"> 2. Continue test with one rack of different columns.

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